




TRANSPORT is the ultimate in impressing convenience. A handheld, disposable impression syringe, TRANSPORT provides more control intraorally compared to bulky impression guns and cartridges, and saves you money on costly mixing tips and wasted impression material.

1



Ensure top slide is open to let air escape (inset). Place TRANSPORT syringe barrels against cartridge and squeeze gun gently. If using a new cartridge, bleed a pea-sized amount of material to ensure even fill (base/catalyst).

2



For 1 prep, fill halfway. For 2 preps, fill the syringe, leaving enough space for the plunger ends to be inserted without premature extrusion at the TRANSPORT mixing section (a little practice is necessary).

Gently insert the plunger to seal the end of the barrels of the TRANSPORT syringe.

3



Slide the mixing tip section all the way across until the tab snaps off.



Depress the plunger to mix and then syringe around the preparation, keeping the tip submerged in the wash material to avoid air entrapment.



**To Pre-Load For Later Use**

Follow steps 1 and 2. Slide mixing tip partway until you hear/feel a click and the O-rings are completely sealed. It is important that you do not snap the tab off during this step. The TRANSPORT syringe is now ready for storage until needed.




TRANSPORT is the ultimate in impressing convenience. A handheld, disposable impression syringe, TRANSPORT provides more control intraorally compared to bulky impression guns and cartridges, and saves you money on costly mixing tips and wasted impression material.

1



Ensure top slide is open to let air escape (inset). Place TRANSPORT syringe barrels against cartridge and squeeze gun gently. If using a new cartridge, bleed a pea-sized amount of material to ensure even fill (base/catalyst).

2



For 1 prep, fill halfway. For 2 preps, fill the syringe, leaving enough space for the plunger ends to be inserted without premature extrusion at the TRANSPORT mixing section (a little practice is necessary).

Gently insert the plunger to seal the end of the barrels of the TRANSPORT syringe.

3



Slide the mixing tip section all the way across until the tab snaps off.



Depress the plunger to mix and then syringe around the preparation, keeping the tip submerged in the wash material to avoid air entrapment.




**To Pre-Load For Later Use**

Follow steps 1 and 2. Slide mixing tip partway until you hear/feel a click and the O-rings are completely sealed. It is important that you do not snap the tab off during this step. The TRANSPORT syringe is now ready for storage until needed.

CE



 167 Central Avenue  
London, ON, Canada, N6A 1M6

For U.S. Distribution:  
Brookfield, CT, USA, 06804

**EC REP** EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands


CLINICIAN'S  
CHOICE®  
DENTAL PRODUCTS INC.

**1-800-265-3444**  
(519) 641-3066  
[www.clinicianschoice.com](http://www.clinicianschoice.com)

092017

CE



 167 Central Avenue  
London, ON, Canada, N6A 1M6

For U.S. Distribution:  
Brookfield, CT, USA, 06804

**EC REP** EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

CLINICIAN'S  
CHOICE®  
DENTAL PRODUCTS INC.

**1-800-265-3444**  
(519) 641-3066  
[www.clinicianschoice.com](http://www.clinicianschoice.com)

092017