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# Evanescce™

Nano-Enhanced Universal Restorative

 First Scientific Dental Materials GmbH  
25335 Elmshorn  
Germany

Distributed by:  
Clinician's Choice Dental Products, Inc  
519.641.3066 **clinicianschoice.com**

 0482

## Instructions for Use

Evanesc™ Nano-Enhanced Universal Restorative is a **visible light-cure nano-filled composite** for anterior and posterior restorations. Beside its low shrink and low abrasion Evanesc™ shows improved handling features and excellent mechanical properties. Evanesc™ is based on urethanedimethacrylate resin and inorganic filler particles <1.0 µm. The total filler content is 82% (wt) and 74% (vol). Radiopacity of Evanesc™ is 200% aluminium (radiopacity of enamel: ~200% aluminium, of dentin: ~100% aluminium). The composite material meets the requirements of **DIN EN ISO 4049**, type 1, class 2, group 1.

### Indications

Suitable for all cavities.

### Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

### Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

### Preparation of Tooth and Color Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate Evanesc™ shade.

Evanesc™ Nano-Enhanced Universal Restorative is shaded according to VITA<sup>®</sup>-Shades. It is also available in special shades and bleach shades.

### Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. For pulp protection, areas close to the pulp should be covered with a thin layer of pulp capping material (e.g. calcium hydroxide liner).

Condition and apply a bonding agent according to the manufacturer instructions. The total-etch technique with an etching gel and subsequent application of a light-cure bonding agent is recommended. Otherwise use a self-etching bonding agent before application of Evanesc™.

## Application and Curing of Evanesc Nano-Enhanced Universal Restorative

Apply Evanesc™ in the selected shade. Evanesc™ Nano-Enhanced Universal Restorative is easy to model. Place it by using a suitable instrument. Matrix strips may be used. In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light-cure resin modified cement or a light-cure flowable composite, followed by the application of Evanesc™.

In case of small cavities Evanesc™ can be applied directly.

For an optimum result apply Evanesc™ in **layers of max. 2mm**.

Light-cure each increment separately.

By using a polymerization unit (wavelength range 400-500 nm) with a **light intensity of at least 1000mW/cm²**, cure each increment as follows:

Lighter shades (e.g. A1, A2, C2) **20 seconds**

Darker shades (e.g. A3.5, B3) **30 seconds**

Hold the light emission window as close as possible to the filling material.

### Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

### Storage

Do not store above 25°C. Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

Single Dose units are for single use only due to hygienic reasons.

### Additional Notes/Warnings

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA<sup>®</sup> is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

### Warranty

First Scientific Dental Materials warrants this product will be free from defects in material and manufacture. First Scientific Dental Materials makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and First Scientific Dental Materials' sole obligation shall be repair or replacement of the First Scientific Dental Materials product.

### Limitation of Liability

Except where prohibited by law, First Scientific Dental Materials will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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