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Template®

Ultra Quick Matrix Material

 First Scientific Dental Materials GmbH
25335 Elmshorn
Germany

Distributed by:
Clinician's Choice Dental Products, Inc
519.641.3066 [clinicianschoice.com](http://www.clinicianschoice.com)

Instructions for Use

Template® Ultra Quick Matrix Material is an extremely fast setting **syringeable one phase impression material** especially for partial impressions of single teeth, groups of teeth or jaw sections and for impressions using the copper band technique.

Template is a medium viscosity, thixotropic and hydrophilic precision impression material based on addition curing vinyl silicones delivered in **automix cartridges**. It can be extruded directly onto the teeth or the tray. Template shows good hydrophilic properties. Thus improves the wettability and yields a better reproduction of subgingival areas.

Template is thixotropic. That means it stays in position in the tray and under pressure it flows. The material meets the requirements of **ISO 4823, type 2** (medium-bodied consistency).

Gun load and extruding

Load the gun with the automix cartridge and mixing cannula. When extruding directly onto teeth use "intraoral tips".

Impression taking

Load the automixed material onto the tray or directly onto the teeth. Insert the tray not later than **30 seconds** (working time at 23°C) after extruding Template into the mouth. Maintain position for at least **50 seconds**.

Note:

Store used cartridge with fixed used mixing cannula. **Do not re-use the cartridge cap!** Closing the cartridge by turning on the cap between use could cause cross contamination of components and thus premature setting.

CASTING

The impression should be stored at normal room temperature (max. 25°C). The impression can be poured **30 minutes** after removal from the mouth. No other time limitation exists.

Electroplating

All impressions can be copper or silverplated.

Disinfection

Impressions can be disinfected with standard disinfection solutions.

Storage

Do not store above 25°C. Do not use after expiry date.

Additional notes:

Do not use with standard (condensation reaction type) silicone impression materials.

Do not use latex gloves, they may interfere with the setting reaction of the impression material.

Warranty

First Scientific Dental Materials warrants this product will be free from defects in material and manufacture. First Scientific Dental Materials makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and First Scientific Dental Materials' sole obligation shall be repair or replacement of the First Scientific Dental Materials product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental Materials will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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