Evanescence® Bulk Cure™

Dual-Cure Bulk Fill Restorative

INSTRUCTIONS FOR USE

Evanescence® Bulk Cure™ is a tooth shaded, radiopaque, dual-cured (auto-cure with light-cure acceleration) composite that offers an unlimited depth of cure in a single application without multiple layers and light-curing when allowed to self-cure. Evanescence Bulk Cure has the handling properties of a low viscosity or “flowable” composite that allows for adaptation to the prepared cavity walls. Evanescence Bulk Cure is also ideal for core buildups, and any application where light transmission is inadequate. It is provided in multiple shades.

INDICATIONS

Direct restoration of Class I, II, III and V cavities (ideally Class I and II), base or liner, core build-ups.

SINGLE-USE MIXING TIPS

Mixing tips are single-use and provided non-sterile. Mixing tips may not function as intended if re-sterilized and may result in an improper procedure and lead to improper function or failure of the device. The mixing tips must be disposed after use on one patient, to avoid cross contamination.

MULTI-USE SYRINGES

Syringes are multi-use devices and are provided non-sterile. Always use a shield, a single use/disposable barrier sheet to prevent exposure or transfer of blood, tissue, or saliva that may cause cross-contamination. Use in accordance with instructions below.

WARNINGS AND PRECAUTIONS

- Keep away from sunlight or ambient light; Replace cap after each use.
- Evanescence Bulk Cure is designed for room temperature use.
- Do not store or expose Evanescence Bulk Cure to elevated ambient temperature or intense light.
- Do not use the composite material in contact with eugenol-containing products, nor store near eugenol-containing materials. Eugenol can impair the polymerization of Evanescence Bulk Cure.

- Evanescence Bulk Cure contains methacrylate/acrylate monomers. Unpolymerized monomers may cause skin sensitization in susceptible persons. Do not use in anyone who has an allergy to acrylics or who has an allergic reaction to the material of this product. In case of contact with skin, wash thoroughly with soap and water.
- For light-cure, use a LED curing light with a minimum output of 1000 mW/cm2. Other light sources or intensities require an adjustment of the cure time. See curing light manufacturer’s instructions.
- Use a single use/disposable barrier sleeve on the syringe to avoid cross-contamination.
- Shelf life of Evanescence Bulk Cure is stated on the syringe labels.
- This product is intended for use by dental professionals only.

COMPOSITIONS

Polymer-based dental restorative material.
Inorganic fillers: average particle size ~0.7µm, ~66% by weight or ~48% by volume.

MIXING WITH MIXING TIPS

Working time at room temperature (22 ± 2°C): > 1 min (without curing light).
Setting time at room temperature (22 ± 2°C): < 10 min (without curing light).

ISO 4049 CLASSIFICATION

Type 1 Class 3.

RADIO-OPAQUE

Evanescence Bulk Cure is radio-opaque.

PROCEDURES FOR USE

Tooth Preparation

- Prepare tooth cavities as needed.
- Apply, and cure a dental bonding adhesive according to its manufacturer’s instruction.

Cavity Restorations

- Twist and remove the cap on the syringe.
- Plug the plunger into the end of the syringe. Push the plunger with smooth and even pressure until pastes move out of the syringe openings. Bleed pastes from two sides on a clean and flat surface to make sure both sides are leveled.
- Manually bend the metal cannula on the mixing tip to the desired angle.
Note: Use a single use/disposable barrier sleeve on the syringe to avoid cross-contamination.
Note: The mixing tips are single use and must be disposed after use on one patient, to avoid cross-contamination.

SAFETY DATA SHEET

ISSUING DATE: 01-Dec-2023
REVISION DATE: 01-Dec-2023
REVISION NUMBER: 1

SECTION 1: IDENTIFICATION

PRODUCT IDENTIFIER

Product Name: Evanescence® Bulk Cure™

OTHER MEANS OF IDENTIFICATION

Product Code(s): 031386100

Synonyms: Bulk Fill

RECOMMENDED USE OF THE CHEMICAL AND RESTRICTIONS ON USE

Recommended use: Dental professional

Restrictions on use: Rx only

DETAILS OF THE SUPPLIER OF THE SAFETY DATA SHEET

Supplier Address/Manufacturer Address:

Distributed by Clinician’s Choice, a division of Den-Mat Holdings, LLC

1017 W. Central Ave., Lompoc, CA 93436 – 800.433.6628

Email: info@clinicianschoice.com

24 Hour Emergency Phone Number: 800.433.6628 M-F

Emergency Telephone: No information available

SECTION 2: HAZARD(S) IDENTIFICATION

CLASSIFICATION: This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

HAZARDS NOT OTHERWISE CLASSIFIED (HNOC): Not applicable.

LABEL ELEMENTS

Hazard statements: Not classified.

OTHER INFORMATION: May be harmful in contact with skin. Causes mild skin irritation.

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

SUBSTANCE: Not applicable.

MIXTURE:

Chemical name	CAS No.	Weight-%	Trade secret
Glass filler	65997-17-3	45-80	*
Triethylene glycol dimethacrylate	109-16-0	2-10	*

*The exact percentage (concentration) of composition has been withheld as a trade secret.

SECTION 12: ECOLOGICAL INFORMATION

Ecotoxicity

Chemical name	Algae/aquatic plants	Fish	Toxicity to microorganisms	Crustacea
Triethylene glycol dimethacrylate	-	LC50: =16.4mg/L (96h, Danio rerio)	-	-

Persistence and degradability: No information available.

Bioaccumulation

Component Information

Mobility: No information available.

Other adverse effects: No information available.

SECTION 13: DISPOSAL CONSIDERATIONS

DISPOSAL METHODS

Waste from residues/unused products: Dispose of in accordance with local regulations.

Dispose of waste in accordance with environmental legislation.

Contaminated packaging: Do not reuse empty containers.

California Hazardous Waste Status: This product contains one or more substances that are listed with the State of California as a hazardous waste.

SECTION 14: TRANSPORT INFORMATION

DOT: Not regulated

SECTION 15: REGULATORY INFORMATION

INTERNATIONAL INVENTORIES

Contact supplier for inventory compliance status

*Contact supplier for details. One or more substances in this product are either not listed on the US TSCA inventory, listed on the confidential US TSCA inventory or are otherwise exempted from inventory listing requirements

US FEDERAL REGULATIONS

SARA 313 Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

SARA 311/312 Hazard Categories Should this product meet EPCRA 311/312 Tier reporting criteria at 40 CFR 370, refer to Section 2 of this SDS for appropriate classifications.

CWA (Clean Water Act) This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

- Attach a mixing tip by securely twisting it onto the head of the syringe, push the plunger with smooth pressure. Do not use the first small amount of material to ensure even mixing.
- Immediately place the mixing tip at the deepest part of the prepared cavity and fill the cavity while slowly moving the tip upward. Always keep the mixing tip in the material to avoid trapping air. Using the metal cannula tip or any other suitable instruments to obtain desired occlusal anatomy as needed.
Note: If a universal composite is desired to be the top layer, only fill the material to the desired portion.
Note: During filling, avoid strong direct light on the mixing chamber as the thickening process may be accelerated.
- Allow the material to self-cure for a minimum of 60 seconds, until it hardens.
- If a universal capping composite is desired, place and cure the universal composite, according to the manufacturer’s instructions.
- Optionally light-cure the occlusal surface for 10 seconds.
- Use an abrasive device to reach the designed occlusal anatomy and polish the surface.
- After use remove the mixing tip and replace the cap securely.

Core Build-Up

- Twist and remove the cap on the syringe.
- Place the plunger into the end of the syringe. Push the plunger with smooth and even pressure until pastes move out of the syringe openings. Remove excess pastes from two sides on a clean and flat surface to make sure both sides are leveled.
- Manually bend the metal cannula on the mixing tip to the desired angle.
Note: Use a single use/disposable barrier sleeve on the syringe to avoid cross-contamination.
Note: The mixing tips are single use and must be disposed after use on one patient, to avoid cross-contamination.
- Attach a mixing tip by securely twisting it onto the head of the syringe, push the plunger with smooth pressure and bleed a small amount of material to ensure even mixing.
- Immediately place the mixing tip at the deepest part of the prepared cavity and fill the cavity while slowly moving the tip upward. Always keep the mixing tip in the material to avoid trapping air. Using the metal cannula tip or any other suitable instruments to obtain desired anatomy as needed.
Note: During filling, avoid strong direct light on the mixing chamber as the thickening process may be accelerated.
- Allow the material to self-cure for at least 60 seconds, until it hardens.

SECTION 4: FIRST-AID MEASURES

DESCRIPTION OF FIRST AID MEASURES

Inhalation: Remove to fresh air.

Eye Contact: Rinse thoroughly with plenty of water, also under the eyelids.

Skin Contact: Wash skin with soap and water.

Ingestion: Rinse mouth.

MOST IMPORTANT SYMPTOMS AND EFFECTS, BOTH ACUTE AND DELAYED

Symptoms: Prolonged contact may cause redness and irritation.

Effects of Exposure: No information available.

INDICATION OF ANY IMMEDIATE MEDICAL ATTENTION AND SPECIAL TREATMENT NEEDED

Note to Physicians: Treat symptomatically.

SECTION 5: FIRE-FIGHTING MEASURES

Suitable Extinguishing Media: Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media: No information available.

Specific hazards arising from the chemical: No information available.

Explosion data:

Sensitivity to mechanical impact: None.

Sensitivity to static discharge: None.

Special protective equipment and precautions for firefighters: Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear.

Use personal protection equipment.

SECTION 6: ACCIDENTAL RELEASE MEASURES

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES

Personal precautions: Ensure adequate ventilation.

METHODS AND MATERIAL FOR CONTAINMENT AND CLEANING UP

Methods for containment: Prevent further leakage or spillage if safe to do so.

Methods for cleaning up: Wash thoroughly after handling.

Prevention of secondary hazards: Clean contaminated objects and areas thoroughly observing environmental regulations.

SECTION 7: HANDLING AND STORAGE

PRECAUTIONS FOR SAFE HANDLING

Advice on safe handling: In syringes, need mixing before use.

CONDITIONS FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES

Storage Conditions: Store in original container. Keep away from incompatible materials.

CERCLA This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302).

US STATE REGULATIONS

California Proposition 65 This product contains the following Proposition 65 chemicals:

Chemical name	California Proposition 65
Titanium dioxide - 13463-67-7	Carcinogen

U.S. STATE RIGHT-TO-KNOW REGULATIONS

Chemical name	N.J.	Mass.	Pa.
Cumene hydroperoxide 80-15-9	X	X	X
Titanium dioxide 13463-67-7	X	X	X
Phenol, 2,6-bis(1,1-dimethylethyl)-4-methyl- 128-37-0	X	X	X
Copper, bis(2,4-pentanedionato-0,0)-, (SP-4-1)- 13395-16-9	X	-	X
Iron oxide 1309-37-1	X	X	X

U.S. EPA LABEL INFORMATION

EPA Pesticide Registration Number: Not applicable

SECTION 16: OTHER INFORMATION

NFPA Health hazards 0 Flammability 0 Instability 0 Special hazards -

HMIS Health hazards 0 Flammability 0 Physical hazards 0 Personal protection X

KEY OR LEGEND TO ABBREVIATIONS AND ACRONYMS USED IN THE SAFETY DATA SHEET

Legend

SVHC: Substances of Very High Concern for Authorization:

PBT: Persistent, Bioaccumulative, and Toxic (PBT) Substances

vPvB: Very Persistent and very Bioaccumulative (vPvB) Substances

STOT: Specific Target Organ Toxicity

ATE: Acute Toxicity Estimate

LC₅₀: 50% Lethal Concentration

LD₅₀: 50% Lethal Dose

- Use an abrasive device, to finish the preparation.
- After use remove the mixing tip and replace the cap securely.

STORAGE

Store at room temperature (36°F/2°C – 77°F/25°C). Do not freeze. Avoid direct sunlight. Avoid excessive heat or light. Use at room temperature.

DISPOSAL

Dispose of used devices which pose a risk of infection according to facility clinical waste procedures and applicable local and state regulations. To dispose of unused material, replace cap and dispose of in accordance with local and state regulations.

NOTICE TO USER IN THE EUROPEAN UNION

Any serious incident that has occurred in relation to the device(s) in which this Instructions for Use applies should be reported to the manufacturer identified in this Instructions for Use and the competent authority of the Member State in which the user and/or patient is established.

Clinician's Choice®

Distributed by Clinician’s Choice, a division of Den-Mat Holdings, LLC
1017 W. Central Avenue, Lompoc, CA 93436 USA
800.433.6628 clinicianschoice.com

753212900 R:0 09/24SN

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

CONTROL PARAMETERS

The following ingredients are the only ingredients of the product above the cut-off level (or level that contributes to the hazard classification of the mixture) which have an exposure limit applicable in the region for which this safety data sheet is intended or other recommended limit. At this time, the other relevant constituents have no known exposure limits from the sources listed here.

Chemical name	ACGIH TLV	OSHA PEL	NIOSH
Glass filler 65997-17-3	TWA: 1 fiber/cm3 respirable fibers: length >5 µm, aspect ratio >=3:1, as determined by the membrane filter method at 400-450X magnification [4-mm objective], using phase-contrast illumination TWA: 5 mg/m3 inhalable particulate matter	-	-

APPROPRIATE ENGINEERING CONTROLS

Engineering controls: Showers, eyewash stations, ventilation systems.

INDIVIDUAL PROTECTION MEASURES, SUCH AS PERSONAL PROTECTIVE EQUIPMENT

Eye/face protection: Wear safety glasses with side shields (or goggles).

Hand protection: Wear suitable gloves.

Skin and body protection: Lab coat. Wear suitable protective clothing.

Respiratory protection: No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

General hygiene considerations: Handle in accordance with good industrial hygiene and safety practice.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

INFORMATION ON BASIC PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Opaque viscous paste

Physical state: Liquid

Color: White to off-white

Odor: No data available

Odor threshold: No data available

PROPERTY / VALUES / REMARKS • METHOD

pH: No data available

pH (as aqueous solution): No data available

Melting point / freezing point: No data available

Initial boiling point and boiling range: No data available

Flash point: > 93 °C / > 199.4 °F

Legend Section 8: Exposure controls/personal protection

TWA: TWA (time-weighted average)

STEL: STEL (Short Term Exposure Limit)

Ceiling: Maximum limit value

Sk*: Skin designation

+ Sensitizers

KEY LITERATURE REFERENCES AND SOURCES FOR DATA USED TO COMPILE THE SDS

Agency for Toxic Substances and Disease Registry (ATSDR)

U.S. Environmental Protection Agency ChemView Database

European Food Safety Authority (EFSA)

EPA (Environmental Protection Agency)

Acute Exposure Guideline Level(s) (AEGLS)

U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act

U.S. Environmental Protection Agency High Production Volume Chemicals

Food Research Journal

Hazardous Substance Database

International Uniform Chemical Information Database (IUCLID)

Japan GHS Classification

Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)

NIOSH (National Institute for Occupational Safety and Health)

National Library of Medicine’s ChemID Plus (NLM CIP)

National Library of Medicine’s PubMed database (NLM PUBMED)

U.S. National Toxicology Program (NTP)

New Zealand’s Chemical Classification and Information Database (CCID)

Organization for Economic Co-operation and Development Environment, Health, and Safety Publications

Organization for Economic Co-operation and Development High Production Volume Chemicals Program

Organization for Economic Co-operation and Development Screening Information Data Set World Health Organization

Symbol Glossary

Symbol	Title	Symbol	Title
	Consult instructions for use		Unique device identifier
	Prescription use only		Batch code
	Temperature limit		Catalogue number
	Keep away from sunlight		
	Do not re-use		
	Use-by date		
	Do not use if package is damaged and consult instructions for use		

SYMPTOMS RELATED TO THE PHYSICAL, CHEMICAL AND TOXICOLOGICAL CHARACTERISTICS

CHARACTERISTICS

Symptoms: Prolonged contact may cause redness and irritation.

ACUTE TOXICITY

Numerical measures of toxicity:

Based on available data, the classification criteria are not met

The following values are calculated based on chapter 3.1 of the GHS document:

ATEmix (oral) 60,833.60 mg/kg

ATEmix (dermal) 2,385.20 mg/kg

COMPONENT INFORMATION

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Triethylene glycol dimethacrylate 109-16-0	= 10837 mg/kg (Rat)	-	-

DELAYED AND IMMEDIATE EFFECTS AS WELL AS CHRONIC EFFECTS FROM SHORT AND LONG-TERM EXPOSURE

Skin corrosion/irritation: Classification based on data available

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