

Precautionary statements

Wear protective gloves and eye protection.

Wear protective clothing.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do so. Continue rinsing.

IF SKIN IRRITATION OCCURS: Get medical advice/attention.

IF ON SKIN: Wash with plenty of water.

Dispose of contents/container in accordance with local/regional/national/international regulations.

SECTION 16: OTHER INFORMATION

The information contained herein is based on our present knowledge and is believed to be correct at the date issued. The information is subject to change without notice. However, all information, statements, recommendations or suggestions shall not constitute any guarantee, representation or warranty, expressed or implied. Final determination whether the product is fit and suitable for a particular purpose and the user's method of use or application is the sole responsibility of the user.

Date of preparation / last revision: 09/08/2016 / 6

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Flam. Liq. 3: Flammable liquids – Category 3

Acute Tox. 4: Acute toxicity – Category 4

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Dam.: 1: Serious eye damage/eye irritation – Category 1

Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A

Skin Sens. 1B: Skin sensation – Category 1B

STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

Asp. Tox. 1: Aspiration hazard – Category 1

***Data compared to the previous version altered.**

Temptation[®] Now

Fast Set Temporary Crown & Bridge Composite

Temptation Now is a fluorescent cold polymerizing paste-paste system for the production of temporary crowns, bridges, inlays and onlays. Temptation Now consists of base paste and catalyst paste.

Temptation Now is available in 3 shades: A1, A2, and A3.

INDICATIONS

Temptation Now is indicated for the fabrication of temporary crowns, bridges, inlays and onlays.

CONTRAINDICATIONS

Temptation Now contains methacrylates, amines, terpenes, benzoyl peroxide, and BHT.

Temptation Now should not be used in cases of known hypersensitivities (allergies) to any of these ingredients.

INSTRUCTIONS

Temptation Now 76 g Cartridges

Insert cartridge into the gun or a compatible dispenser with a suitable piston. Due to technical reasons, slight differences in the filling level between both cartridge chambers may occur. In this case remove cartridge cap and squeeze out the material until it is delivered evenly from both chambers. Then, insert a mixing tip and lock with a 90° clockwise turn. By pressing the lever the material is extruded and automatically mixed in the correct ratio. The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.

Do not remove the mixing tip from the cartridge after use. Replace the used mixing tip with a new one just before the next application. Check the cartridge opening to ensure material flowability when changing the mixing tips.

Temptation Now 8 g Syringes

Due to technical reasons slight differences in the filling level between the syringe chambers can occur prior to the first use. Remove the cap from the syringe after turning it counterclockwise and squeeze the material out until the material is delivered evenly from both chambers. Attach a mixing tip.

Match the nooks on the syringe to the mixing tip and lock the mixing tip into place with a 90° clockwise turn. The material is automatically mixed in the correct ratio when dispensing and can be applied directly. All syringes come with a detachable piston. The piston can simply be pulled from the syringe and replaced without causing any material leakage. The plugs in the syringe chambers will stay in place preventing any cross-contamination of the material in the syringe.

APPLICATION

Level out any undercuts in the preparation. Use only eugenol-free cements. In cases where no pre-fabricated crowns or similar restorations are used, prepare an impression by levelling undercuts and cutting drainage channels. If the wall of the temporary restoration is too thin, enlarge the impression, especially removing septums in the interdental embrasures. Then clean and air dry the impression. Apply Temptation Now directly into the impression. Do not overfill.

Temptation Now attains an elastic phase after 1-1.5 min from beginning of mixing (approx. 0.5-1 min after inserting into the mouth).

Remove the temporary from the mouth no later than 1.5 min (elastic phase) after mixing, to ensure a trouble-free removal.

The setting process can additionally be monitored by checking the excess material in the mouth.

The setting time of non-used material (in mixing tips) cannot be compared with the setting time of the material in the mouth.

At a room temperature of 23°C and 50% humidity working times are as follows:

Min:s 0:00	Min:s 0:30	Min:s 1:00	Min:s 1:30	Min:s 4:00
Application in the Mouth	Setting in the Mouth	Removal	Complete Setting	Finishing

Slight differences may result from variations in room temperatures. After approx. 2-3 min, remove the temporary from the impression and wipe dry. If necessary, use a solvent, e.g. ethyl alcohol. This removes the thin inhibition layer caused by oxygen and allows easier finishing. The surface of the temporary restoration should appear silky-matt. Carefully trim any excess material and roughen areas near the gingiva. Cement the temporary restoration using a eugenol-free temporary cement (e.g. Cling² Temporary Cement, Clinician's Choice)

CAUTION

- Before making the temporary, the prepared cores should be insulated for desensitization (fluoride varnishes, core varnishes, etc.).
- In order to avoid fractures or injuries from debris, the temporary should not be exposed to excessive loads (molar area).
- Linings and cores built-up with resin-based materials (composites) should be sufficiently insulated (e.g. glycerine) for easy separation of the temporary from composite resin material.
- Inadequate hygiene or intensive use of some mouthwashes may lead to discoloration.
- Finishing should only be started after complete setting of the material (approx. 4 min after mixing).

REPAIRS

If the temporary breaks despite Temptation Now's extreme stability, the following procedures are recommended:

- Air bubbles: fill with Temptation Now
- In case of fracture immediately after preparing the temporary: re-connect the fracture with freshly mixed Temptation Now
- Restoration of already used temporaries: Roughen the fracture with a bur and place undercuts. Apply a bonding agent to the prepared fractures and cure according to the respective instructions for use. Apply a flowable composite to the prepared areas. Bond the fractured pieces together and cure according to the instructions for use.

STORAGE

- Store at temperatures between 39°F-73°F (4°C-23°C). Do not use after expiry date.
- Do not expose Temptation Now to direct sunlight or the chairlight.
- This material has been developed solely for use in dentistry.
- Processing should be done strictly according to the instructions for use.

Clinician's Choice Dental Products, Inc. recognizes its responsibility to replace products if proven to be defective. Clinician's Choice does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

Keep this material out of reach of children.

For dental use only.



Distributed by Clinician's Choice, a division of Den-Mat Holdings, LLC
1017 W. Central Avenue, Lompoc, CA 93436 USA
800.433.6628 [clinicianschoice.com](http://www.clinicianschoice.com)



Temptation[®] Now

Fast Set Temporary Crown & Bridge Composite

800.433.6628
www.clinicianschoice.com

SAFETY DATA SHEET

SECTION 1: IDENTIFICATION

PRODUCT IDENTIFIER

Trade name: Temptation Now Fast Set Temporary Crown & Bridge Composite

Application of the substance / the mixture: Temporary Crown & Bridge Material

DETAILS OF THE SUPPLIER OF THE SAFETY DATA SHEET

Manufacturer/Supplier:

Distributed by Clinician's Choice, a division of Den-Mat Holdings, LLC

1017 W. Central Avenue, Lompoc, CA 93436 USA

For emergencies or product information 800.433.6628

or email info@clinicianschoice.com

SECTION 2: HAZARDS IDENTIFICATION

Health effects are based on the unmixed/uncured material.

Under FDA regulations this product is classified as a medical device.

Medical devices are exempt from the labeling provisions of GHS according to the Hazard Communication Standard (HCS).

Classification of the substance or mixture: The product is classified according to the Globally Harmonized System (GHS).



Skin Sens. 1B H317 May cause an allergic skin reaction.

Label elements:

Hazard pictograms: GHS07

Signal word: Warning

Hazard statements: May cause an allergic skin reaction.

Precautionary statements:

Wear protective gloves and eye protection.

Wear protective clothing.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

IF SKIN IRRITATION OCCURS: Get medical advice/attention.

IF ON SKIN: Wash with plenty of water.

Dispose of contents/container in accordance with local/regional/national/international regulations.

Other hazards: Use only for intended scope of application

Results of PBT and vPvB assessment:

PBT: Not applicable.

vPvB: Not applicable.

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

Chemical characterization:

Mixtures

Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Polyesterdimethacrylate	Skin Irrit. 2, H315; Eye Irrit. 2A, H319	10-25%
Catalyst	Eye Dam. 1, H318 Acute Tox. 4, H302 Flam. Liq. 3, H226	2.5-5%
Co-catalyst	Asp. Tox. 1, H304 Skin Irrit. 2, H315; Eye Irrit. 2A, H319; STOT SE3, H335	1-2.5%

SECTION 4: FIRST AID MEASURES

DESCRIPTION OF FIRST AID MEASURES

After inhalation: Seek medical treatment in case of complaints.

After skin contact: Immediately wash with water and soap and rinse thoroughly.

After eye contact: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

After swallowing: If symptoms persist consult doctor.

INFORMATION FOR DOCTOR

Most important symptoms and effects, both acute and delayed: No further relevant information available.

Indication of any immediate medical attention and special treatment needed: No further relevant information available.

SECTION 5: FIRE FIGHTING MEASURES

EXTINGUISHING MEDIA

Suitable extinguishing agents: CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

Special hazards arising from the substance or mixture: No further relevant information available.

ADVICE FOR FIREFIGHTERS

Protective equipment: No special measures required.

SECTION 6: ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures: Not required.

Environmental precautions: No special measures required.

Methods and material for containment and cleaning up: Pick up mechanically.

REFERENCE TO OTHER SECTIONS

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7: HANDLING AND STORAGE

HANDLING

Precautions for safe handling: For dental use only.

Information about protection against explosions and fires: No special measures required.

Conditions for safe storage, including any incompatibilities

STORAGE

Requirements to be met by storerooms and receptacles: See storage instructions on package

Information about storage in one common storage facility: Not required

Further information about storage conditions: Keep receptacle tightly sealed.

Specific end use(s) No further relevant information available.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

Additional information about design of technical systems: No further data; see item 7.

CONTROL PARAMETERS

Components with limit values that require monitoring at the workplace: The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

Additional information: The lists that were valid during the creation were used as basis.

EXPOSURE CONTROLS

Personal protective equipment:

General protective and hygienic measures: The usual precautionary measures for handling chemicals should be followed.

Breathing equipment: Not required.

Protection of hands: Protective gloves



Material of gloves: The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and; therefore, has to be checked prior to the application.

Eye protection: Safety glasses

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

GENERAL INFORMATION

Appearance:

FORM: Pasty

COLOR: Light brown

Odor: Characteristic

Odor threshold: Not determined.

pH-value: Not determined.

CHANGE IN CONDITION

Melting point/Melting range: Undetermined.

Boiling point/ Boiling range: Undetermined.

Flash point: Not applicable.

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Density: Not determined.

Relative density Not determined.

Vapor density

Not applicable.

Evaporation rate

Not applicable.

Solubility in / Miscibility with Water: Not miscible or difficult to mix.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

Not determined.

Dynamic:

Not determined.

Kinematic:

Not determined.

Other information

No further relevant information available.

SECTION 10: STABILITY AND REACTIVITY

REACTIVITY: No further relevant information available.

CHEMICAL STABILITY

Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.

Possibility of hazardous reactions: No dangerous reactions known.

Conditions to avoid: No further relevant information available.

Incompatible materials: No further relevant information available.

Hazardous decomposition products: No dangerous decomposition products known.

SECTION 11: TOXICOLOGICAL INFORMATION

INFORMATION ON TOXICOLOGICAL EFFECTS

Acute toxicity:

Primary irritant effect:

ON THE SKIN: Testing according to ISO 10993-1/10993-10 showed no skin irritation for this product when fully cured

ON THE EYE: No irritating effect.

Sensitization: Sensitization possible through skin contact.

Additional toxicological information: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

CARCINOGENIC CATEGORIES

IARC (International Agency for Research on Cancer):

fumed silica

128-37-0 2,6-di-tert-butyl-p-cresol

NTP (National Toxicology Program): None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration): None of the ingredients is listed.

SECTION 12: ECOLOGICAL INFORMATION

TOXICITY

Aquatic toxicity: No further relevant information available.

Persistence and degradability: No further relevant information available.

BEHAVIOR IN ENVIRONMENTAL SYSTEMS

Bioaccumulative potential: No further relevant information available.

Mobility in soil: No further relevant information available.

RESULTS OF PBT AND vPvB ASSESSMENT

PBT: Not applicable.

vPvB: Not applicable.

Other adverse effects: No further relevant information available.

SECTION 13: DISPOSAL CONSIDERATIONS

WASTE TREATMENT METHODS

Recommendation: Smaller quantities can be disposed of with household waste.

Uncleaned packagings:

Recommendation: Disposal must be made according to official regulations.

SECTION 14: TRANSPORT INFORMATION

UN-NUMBER

DOT, ADR, AND, IMDG, IATA Void

UN PROPER SHIPPING NAME

DOT, ADR, AND, IMDG, IATA Void

TRANSPORT HAZARD CLASS(ES)

DOT, ADR, AND, IMDG, IATA Void

Class

DOT, ADR, AND, IMDG, IATA Void

PACKING GROUP

DOT, ADR, AND, IMDG, IATA Void

ENVIRONMENTAL HAZARDS

Marine pollutant: No

Special precautions for user Not applicable.

TRANSPORT IN BULK ACCORDING TO ANNEX II

Of MARPOL73/78 and the IBC Code Not applicable.

UN "Model Regulation": Void

SECTION 15: REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

SARA

Section 355 (extremely hazardous substances): None of the ingredients are listed.

Section 313 (Specific toxic chemical listings): None of the ingredients are listed.

TSCA (Toxic Substances Control Act): Glas Powder

PROPOSITION 65

Chemicals known to cause cancer: None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for females: None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for male: None of the ingredients are listed.

Chemicals known to cause development toxicity: None of the ingredients are listed.

CARCINOGENICITY CATEGORIES

EPA (Environmental Protection Agency): None of the ingredients are listed.

TLV (Threshold Limit Value established by ACGIH): Contact Clinician's Choice for more information.

NIOSH-Ca (National Institute for Occupational Safety and Health): None of the ingredients are listed

GHS label elements: The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms: GHS07

Signal word: Warning

Hazard statements: May cause an allergic skin reaction.

CONTINUED...