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Affinity™

Hydroactive Impression Material

 First Scientific Dental Materials GmbH
25335 Elmshorn
Germany

Distributed by Clinician's Choice, a division of Den-Mat Holdings, LLC
1017 W. Central Avenue, Lompoc, CA 93436 USA
800.433.6628 clinicianschoice.com

Instructions for Use

INTENDED USE: A vinyl poly siloxane material to be used as a medical device for the purposes of creating, via impression, a reproduction of tooth and gum structure.

Hydroactive Affinity™ V.P.S. is the first 3rd generation vinyl poly siloxane impression material. Affinity's advanced 3rd generation chemistry provides true hydroactive properties, higher tear strength and dimensional stability. In addition, each viscosity is matched to all others resulting in a superior impression, regardless of the technique utilized. Affinity also exhibits the most efficient working and set times to minimize chair time and the potential for inaccuracies. Affinity's dispensing is even advanced in that Affinity flows more easily through a conventional mixing tip as compared to previous generation impression materials, regardless of the viscosity. 3rd generation Affinity will consistently provide you with the most efficient and accurate impression.

Instructions for Automix Cartridges

- 1) Insert the flange at the rear of the cartridge into the space provided at the front of the impression gun – depress the gun lever to lock the cartridge in place.
- 2) Twist off and dispose of the cartridge shipping cap.
- 3) Attach a mixing tip and twist 1/4 turn to lock it into position; the mixing is done automatically inside the mixing tip.
- 4) Before filling the tray, it is recommended to discard the initially mixed (pea-sized) amount of the material. Check to ensure both base and catalyst are expressing freely.
- 5) Gently squeeze the impression gun handle to mix and dispense the impression material. Release handle to stop the flow.
- 6) Leave the used automix tip on the cartridge as a self sealer after each use. For subsequent dispensing remove and discard the sealer tip. Proceed with steps 3 through 6.

Mixing/Setting Times

AFFINITY V.P.S. HYDROACTIVE	COLOR	AVAILABLE WORKING TIME	MINIMUM SET TIME	INTRAORAL TIME	MAXIMUM TOTAL CURE TIME FROM START OF MIX
LIGHT H.F. REGULAR SET	LIGHT BLUE	1:45 min	2:30 min	4:15 min.	
<small>Meets the requirements of ISO 4823 (type 3).</small>					
LIGHT H.F. FAST SET	LIGHT BLUE	1:10 min	1:45 min	2:55 min.	
<small>Meets the requirements of ISO 4823 (type 3).</small>					
HEAVY BODY REGULAR SET	TEAL	1:45 min	2:30 min	4:15 min.	
<small>Meets the requirements of ISO 4823 (type 2).</small>					
HEAVY BODY FAST SET	DK. GREEN	1:10 min	1:45 min	2:55 min.	
<small>Meets the requirements of ISO 4823 (type 2).</small>					
INFLEX REGULAR SET	DK. BLUE	1:45 min	2:30 min	4:15 min.	
<small>Meets the requirements of ISO 4823 (type 1).</small>					
INFLEX FAST SET	DK. BLUE	1:10 min	1:45 min	2:55 min.	
<small>Meets the requirements of ISO 4823 (type 1).</small>					
QUICK BITE	PURPLE	0:15 min	0:45 min	1:00 min.	

AFFINITY FAST SET MATERIALS ARE RECOMMENDED FOR SINGLE PREPARATIONS ONLY.
ALL WORKING TIMES ARE CALCULATED FOR 23°C ROOM TEMPERATURE.

Storage

Do not store above 25°C. Do not use after expiration date.

Warranty

First Scientific Dental Materials warrants this product will be free from defects in material and manufacture. First Scientific Dental Materials makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and First Scientific Dental Materials' sole obligation shall be repair or replacement of the First Scientific Dental Materials product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental Materials will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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