Zementieren

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 Clinician's Choice, a division of Den-Mat Holdings, LLC
1017 W. Central Avenue, Lompoc, CA 93436 USA
800.433.6628 **clinchianschoice.com**

 **0482**

988234 R-2 02/25CM

R_x only

MD

Instructions for Use

Cling2® is an **automix non-eugenol temporary crown and bridge cement** suitable for cementing temporary crowns or bridges as well as for trial cementing of permanent restorations.

It shows excellent adhesion to temporary crown and bridge materials, but it can easily be removed from prepared tooth stumps. Cling2 is delivered in minimix syringes and/or automix cartridges so it is very easy and comfortable to use. Due to its non-eugenol formulation it does not influence the setting of resin-based filling and luting materials. The temporary crown and bridge cement is based on ZnO/organic acids. Especially in cases of allergies to eugenol it is strongly recommended to use this eugenol-free type of cement.

Cling2 meets the requirements of **DIN EN ISO 3107** and **ADA/ANSI Specification No. 30**.

Application

Preparing the 10 mL or 5 mL Minimix-Syringe

Remove the cap of the minimix-syringe and throw it away (**do not use it again!**). It is replaced by a special 1:1 mixing cannula. Turn the cannula 90° until it locks into position. The material is now ready for application.

Note:

Store used syringe with fixed used mixing cannula.

Discard the first 2-3 mm of the extruded material. This must be done for each new mix.

Preparing the 25 mL Automix Cartridge:

Application gun and cartridge

Load the application gun with 1:1 plunger and remove and dispose the closure-cap. **Do not re-use the cap!** Install the mixing cannula onto the cartridge and extrude needed material. After each use leave the used mixing cannula in place to seal the cartridge.

Note:

The initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix.

Cementing

Cling2 can be applied directly into the temporary crown or bridge under slight pressure from the mixing cannula.

Apply a thin layer on the inside of the temporary crown or bridge and seat no later than **1 minute** (working time at 23°C) cautiously but firmly over the prepared tooth structures.

Note:

Do not dry stumps, best moisten them with some water. A humid milieu (humid stumps) additionally facilitates the removal of the temporary crown or bridge later on.

Cling2 sets in the mouth within **2-3 minutes**. After Cling2 has set remove any excess.

Precautions

1. Do not use Cling2 in cases where known allergies to the components exist.
2. Do not use Cling2 if in some cases individual contact allergies may exist.
3. In case of contact with skin and eyes rinse immediately with plenty of water and consult a physician if necessary.
4. If swallowed seek medical advice.

Storage

Store at 8 - 25°C. Do not use after expiry date.

Warranty

First Scientific Dental Materials warrants this product will be free from defects in material and manufacture. First Scientific Dental Materials makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and First Scientific Dental Materials' sole obligation shall be repair or replacement of the First Scientific Dental Materials product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental Materials will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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